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IFU-ANL v01/2019-PT

INSTRUCTIONS FOR THE USE OF REPLICAS

1. Product description

Replicas are machined refractory stainless steel or medical class titanium elements designed to support high temperatures in the event of welding.

2. Planned use:

They are used to replicate the transepithelial/digital implant of plaster laboratory model and so their name.

3. Directions for use

Before any manipulation, make sure that the anti-rotation and retaining parts are firmly connected.

Before tightening, make sure that the replica and the prosthesis coincide in size and connection type.

Replicas are machined with furrows or facets which, once the plaster sets, prevent any movement and fix them securely in place.

4. Cleaning, disinfection and sterilisation

In the event of any packaging defect, the product should not be used and the manufacturer notified immediately.

All BIONER replicas are supplied clean and appropriately packaged. As they are not inserted in the mouth and their use is restricted to the laboratory, they require no sterilisation.

5. Warnings

Do not reuse any replica, they are products for single use only; reuse could alter their connections and technical specifications.

6. Precautions

It is especially important to ensure that each reference is used in combination with dimensionally compatible elements as indicated on the product label

7. Compatibility


BIONER replica family references must be used in conjunction with compatible components. Carefully examine the components available in the company product catalogue for suitability before insertion.

8. Validity

These instructions for use correspond to the current version of the replica family references indicated at the foot of the document.

9. Storage

BIONER products do not require any special storage conditions even though it is recommended not to expose them to conditions that could affect the product, its packaging or label, such as extreme temperature and humidity, knocks, etc.

| ISO 15223 Symbol | Description |
|---|--|
|  | Indicates the direction for opening the blister |
|  | BIONER SA logotype |
| REF | Catalogue reference |
| LOT AOFAAXXXXX | Batch code |
|  | Date of Manufacture |
|  | Refer to the instructions for use |
|  | Do not use if package is damaged |
| Family | Name of the product family |
| Mat = Stainless steel/Titanium | General indication of the material used |
| q=1 | Indicates the number contained in each blister pack |
|  | Non-sterile |
|  | QR code (reference and batch) |
|  | Do not reuse |
|  | Consult warnings and precautions in instructions for use |
|  | Manufacturer |
| CE | EC Marking |

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