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## INSTRUCTIONS FOR USE OF SLEEVES FOR GUIDED SURGERY

### 1. Product description

The BIONER SA family of guided surgery sleeves are machined from single use biocompatible Titanium. They are cylindrical in shape and available in different inside diameters and lengths to simplify fully guided drilling and implant placement. This product family also includes mounting pins to hold the sleeve in place during surgery.

### 2. Planned use:

The elements included in this family are used to guide computer aided drilling and implant placement processes based on precise preoperative planning. The guided surgery sleeves define the position and three-dimensional orientation of the implants and so contribute to minimising problems during surgery.

### 3. Indications for use

Before use, make sure that the selected elements are suited to the surgical drills to be used and the size of the dental implant to be placed.

### 4. Contraindications

Allergy or hypersensitivity to the chemical components of the materials used.

### 5. Side effects and interactions

Follow the instructions for use recommended by the manufacturer. Incorrect use of the product may cause injury to the patient,

### 6. Cleaning, disinfection and sterilization

The product should not be used and the manufacturer notified immediately in the event of any packaging defect.

Guided surgery sleeves are marketed as non-sterile products. They should be cleaned, disinfected and sterilised by the healthcare professional before use.

Bioner SA always recommends performing previous cleaning with products commonly used in dentistry, with emphasis on ultrasonic baths or automatic cleaning and disinfection methods before proceeding with sterilisation.

Remove the product from its original packaging and place it in a suitable container before sterilisation.

Sterilisation should be performed in sterilisation bags or boxes in a wet steam autoclave at 121°C for 15 minutes or 134°C for at least 3 minutes, as per UNE-EN-ISO 17665-1:2007.

Do not use the product after sterilisation if there is any detectable deterioration of the container.

Guided surgery sleeves, as well as all accessories produced by Bioner SA, should not be reused; they are products for single use only.

## 7. Warnings

Guided surgery sleeves should not be reused; they are products for single use only, as any reuse could alter their technical specifications.

The healthcare professional using the product must have experience in using this product safely and in accordance with instructions for use provided by BIONER SA, the said professional being liable for the suitability of the component to the patient's needs.

## 8. Precautions

It is especially important to ensure that each reference is used in combination with dimensionally compatible elements as indicated on the product label.

## 9. Compatibility


The references of the BIONER family of guided surgery sleeves must be used in conjunction with compatible components. Carefully examine the components available in the company product catalogue for suitability before insertion.

## 10. Validity

These instructions for use correspond to the current version of references of the family of guided surgery sleeves as indicated at the foot of the document

## 11. Storage

BIONER products do not require any special storage conditions even though it is recommended not to expose them to conditions that could affect the product, its packaging or label, such as extreme temperature and humidity, knocks, etc.

ISO 15223 Symbol	Description
	Indicates the direction for opening the blister
 <b>BIONER</b> SISTEMAS IMPLANTOLÓGICOS	BIONER SA logotype
<b>REF</b>	Catalogue reference
<b>LOT</b> AOFAAXXXX	Batch code
 DD/MM/YY	Date of Manufacture
	Refer to the instructions for use

	Do not use if package is damaged
Family	Name of the product family
Mat = Titanium	General indication of the material used
q=x	Indicates the amount contained in each blister
	Non-sterile
	QR code (reference and batch)
	Do not reuse
	Consult warnings and precautions in instructions for use
	Manufacturer
	CE mark

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