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INSTRUCTIONS FOR THE USE OF A KIT FOR DENTAL GRAFTS

1. Product description

The kit for dental grafts and its replacement components consist of:

- Mounting accessories: elements used directly to ensure bonding of the graft during surgery and/or buccal reconstruction. Specifically, these elements are:
 - E- INJ-CHIC. Tack for the fixation of membranes.
 - F- INJ-MDM. DM self-tapping micro screw for the fixation of soft tissues.
- Placement devices: different tools (handles and screwdrivers) designed for placement of mounting accessories.
- Support device: box for storage of different components.

2. Planned use:

The kit for dental grafts and its replacement parts are designed for use in oral surgical procedures that require the fixation of mucous, bone and regeneration membrane grafts, as well as for stabilising the edges in vestibuloplasty. The mounting accessories may also be used as complements for sutures.

The kit for dental grafts should be used under controlled environmental conditions

The fixing elements are indicated for temporary fixing and should be removed once the reconstruction has been consolidated.

3. Contraindications

All mounting accessories are manufactured from medical grade titanium.

All placement devices (tools) are manufactured from medical grade stainless steel, with the exception of the tips of the tacks which are medical grade titanium.

Contraindications: allergy or hypersensitivity to the chemical components of the materials used.

Use of the product is also contraindicated in cases of active infection, in patients showing previous sensitivity to titanium or stainless steel, as well as patients with disorders that could make them incapable of following the physician's instructions. Use is not recommended in the event of systemic, bone or metabolic disorders that could alter the result of the reconstruction, as well as in the case of drug or alcohol abuse or medical conditions under treatment with steroids or which affect blood coagulation.

4. Side effects and interactions

Follow the instructions for use recommended by the manufacturer. Incorrect use of the product may cause injury to the patient.

5. Cleaning, disinfection and sterilisation.

Kits for dental grafts are supplied unsterilised. Before using them, the healthcare professional must perform cleaning, disinfection and sterilisation. Bioner SA always recommends performing previous cleaning of all components with products commonly used in dentistry, with emphasis

on ultrasonic baths or automatic cleaning and disinfection methods before proceeding with sterilisation.

Wear gloves when cleaning the devices and always use non-abrasive disinfectant agents. The preferred sterilisation technique is wet steam autoclave in accordance with general instructions provided by BIONER S.A.

Sterilisation of the entire kit should be performed in sterilisation bags in a wet steam autoclave at 121°C for 15 minutes or 134°C for at least 3 minutes, as per UNE-EN-ISO 17665-1:2007.

Remove the kit from the sterilisation bag immediately before use. Remove the cover from the kit. If necessary, fit the handpiece with a suitable tip, make sure the connection is correct and that the fixing device to be used is firmly held in the placement tool.

Do not use the product after sterilisation if there is any detectable damage.

Multiple re-sterilisation of the mounting elements may alter their working capacity because of successive changes in temperature.

Remove the cover from the organiser to dispense screws and instruments, placing the different items in their corresponding place. Locate the placement elements to be required for the operation, once they are clean and dry. Place the number of mounting points required for the operation with planning for any possible incidents and needs. Replace the cover on the box and place the prepared kit in a sterilisation bag of suitable size. Close the bag and label it with the sterilisation date, patient code and any other routine information.

Fixing accessories should not be reused; they are products for single use only.

6. Warnings

Before use, make sure that the material is in perfect condition and there are no loose parts to prevent the inhalation or swallowing of one or more components. Furthermore, it is highly recommended to make sure that the product is clean and disinfected and the part in contact with the patient is clean and free of any residue.

The healthcare professional using the product must have experience in using this product safely and in accordance with instructions for use provided by BIONER SA, the said professional being liable for the suitability of the component to the patient's needs.

The storage box for the tools and different fixing accessories may be reused, but not the fixing accessories themselves (tacks and micro-screws).

7. Precautions

It is especially important to ensure that each reference is used in combination with dimensionally compatible elements as indicated on the product label.

Inspect the placement elements before and after each procedure to ensure they are in good working condition. Defective instruments must not be used and should be replaced.

8. Compatibility

The references for BIONER SA. kits for dental grafts should be used in conjunction with compatible components, and the size of the accessories selected in accordance with the surgical procedure and type of graft to be performed.

Carefully examine the components available in the company product catalogue for suitability before insertion.

Make sure to use the correct tool for each type of mounting element.

The dimensions ensure compatibility with the components matching each specific reference. The compatibility of each reference is guaranteed by BIONER SA in-process controls.

9. Validity

These instructions for use correspond to the current version of references for the kit for dental grafts family as indicated at the foot of the document.

10. Storage

Elements of this family must be preserved and stored in non-extreme environmental conditions, avoiding prolonged exposure to sunlight.

ISO 15223 Symbol	Description
	Indicates the direction for opening the blister
	BIONER SA logotype
REF	Catalogue reference
LOT AOFAAXXXX	Batch code
	Date of Manufacture
	Refer to the instructions for use
	Do not use if package is damaged
Family	Product family name / additional info
Mat = Ti and Stainless steel	General indication of the material used
q=1	See product description
	Non-sterile
	QR code (reference and batch)
	Do not reuse (except placement and support devices)
	Consult warning and precautions in instructions for use
	Manufacturer
	EC Marking

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