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INSTRUCTIONS FOR USE ABUTMENTS

1. Product description

The abutments are solid elements of titanium or Cr-Co / POM, which can be screwed both directly onto the implant, and onto a transepithelial connector.

2. Intended use

The abutments constitute the framework on which the dental prosthesis rests, and normally, they are adapted in the laboratory before their use in clinic.

The Chromium-Cobalt / POM abutments undergo a transformation of the POM part into dental metal in the prosthetic laboratory before being incorporated into the mouth.

3. Directions for use

The abutments are metallic elements to perform prosthetic restorations on the implants directly or on transepithelial connectors.

It is important to ensure the choice of the abutment regarding the model (compatibility, height and angulation).

There are provisional abutments that are intended to support temporary prostheses or immediate loads. The recommended time of use is 30 days, and after the time prescribed by the doctor they must be replaced by the definitive ones

4. Side effects

Allergy and hypersensitivity to the chemical components of the materials used (medical grade titanium).

5. Side effects and interactions, complications with implants

It is advisable to follow the instructions for use recommended by the manufacturer. Improper use of the product may result in patient damage.

6. Cleaning and disinfection and sterilisation.

The abutments are marketed without sterilisation. Before using them, the professional must perform cleaning, disinfection and proceed to sterilisation.

Bioner SA recommends to always carry out a previous cleaning with products of use in the dental sector, relying on ultrasonic baths or automatic cleaning, disinfection methods and then proceed to sterilisation.

Prior to sterilising, the products in their original packaging must be transferred to another one suitable for said purpose.

Sterilisation must be carried out in a steam autoclave at 121°C for 15 minutes or 134°C at a minimum of 3 minutes, according to UNE-EN-ISO 17665-1: 2007, in sterilisation pouches or trays.

Do not use the product after sterilisation if any deterioration in the package is observed.

The abutments, as well as all the attachments produced by Bioner SA, must not be reused; they are single-use products.

7. Warning

Before use, it must be checked that the material is in perfect condition and there are no disassembled parts, to avoid aspiration or swallowing of part or all components. Furthermore, it is especially recommended to check that the product is clean and disinfected and that the patient connection area is clean and free of debris.

The healthcare professional who uses the reference product must have experience to use this product in a safe and adequate way to the instructions provided by BIONER SA, being the responsibility of the same professional to apply this element to the needs of the patient.

8. Precautions

It must be especially taken into account that each reference must be used in combination with dimensionally compatible elements and which are indicated on the product label itself.

9. Compatibility

BIONER SA abutment family references must be used in conjunction with compatible components. Before fitting, it is necessary to review the components available in the company's product catalogue to make the correct choice


The correct dimensions of the product ensure compatibility with the components with which each specific reference will be combined. The guarantee of the compatibility of each reference is based on the controls of the manufacturing process performed by BIONER SA.

10. Validity

These instructions for use correspond to the current version of the abutment family references, indicated at the bottom of the document.

11. Storage

Elements of this family must be preserved and stored in non-extreme environmental conditions, avoiding prolonged exposure to sunlight.

ISO symbol 15223	Description
	Indicator of the blister opening direction
	BIONER SA logo
REF	Catalogue reference
LOT AOFAAXXXX	Batch Code
 DD/MM/YY	Manufacturing date

	See instructions for use
 Don't use if package is damaged.	Do not use if package is damaged
Family	Product family name / additional info
Mat = Ti (or Cr-Co / POM)	Generic code of the construction material
q=1	Indicative of the amount contained in each blister
	Non-sterile
	QR code (reference and lot)
	Do not re-use
	Consult warning and precautions in instructions for use
	Manufacturer
	EC Marking

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