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IFU-TFR v01/2019-PT

OBSOLETE

INSTRUCTIONS FOR THE USE OF A KIT FOR DRILL STOPS

25. Product description

BIONER SA. kits for drill stops contain machined grade 5 titanium biocompatible, reusable, cylindrical elements of different diameters and lengths, whose inside diameter is suitable for placement on the drill.

They are available on the market in kits that are appropriately identified according to their planned use.

26. Planned use:

The purpose of drill stops is to ensure the drill depth does not exceed the size of the dental implant. They are for optional use by the healthcare professional.

27. Directions for use

Before use, make sure that the selected drill stops correspond to the surgical drills being used and the depth of the dental implant to be placed. They are colour-coded to ensure correct use. The colour indicates the drill diameter and the laser marks on the drill show the drill depth.

28. Contraindications

Allergy or hypersensitivity to the chemical components of the materials used; medical grade Ti.

29. Side effects and interactions

Follow the instructions for use recommended by the manufacturer. Incorrect use of the product may cause injury to the patient.

30. Cleaning, disinfection and sterilisation

Kits for drill stops are delivered unsterilised. The healthcare professional should clean, disinfect and sterilise the devices before use.

Bioner always recommends performing previous cleaning with products commonly used in dentistry, with emphasis on ultrasonic baths or automatic cleaning and disinfection methods before proceeding with sterilisation.

Sterilisation should be performed in sterilisation bags in a wet steam autoclave at 121°C for 15 minutes or 134°C for at least 3 minutes, as per UNE-EN-ISO 17665-1:2007.

Remove the kit from the sterilisation bag immediately before use.

Do not use the product after sterilisation if there is any detectable deterioration of the container.

31. Warnings

Drill stops may be reused after cleaning, disinfection and sterilisation.

The healthcare professional using the product must have experience in using this product safely and in accordance with instructions for use provided by BIONER SA, the said professional being liable for the suitability of the component to the patient's needs.

32. Precautions

It is especially important that each reference be used in combination with dimensionally compatible elements.

33. Compatibility

The compatibility of every reference is guaranteed by the stringent controls BIONER SA applies to all devices. Refer to the drilling guide in the BIONER catalogue.

34. Presentations


Drill stops are available in kits for procedures.

35. Validity

These instructions for use correspond to the current version of references for the kit for drill stops family as indicated at the foot of the document.

36. Storage

Elements of this family must be presented and stored in non-extreme environmental conditions, avoiding prolonged exposure to sunlight.

ISO 15223 Symbol	Description
	BIONER SA logotype
REF	Catalogue reference
LOT AOFXXXXX	Batch code
 DD/MM/YY	Date of Manufacture
	Refer to the instructions for use
 Don't use if package is damaged.	Do not use if package is damaged
Name of the KIT	Name identifying the KIT
Components	Description of each of the components on the box.
	Non-sterile
Sterilisation method	Sterilisation method, temperature and time
	Precautions
	Manufacturer

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Name BIONER, S.A.