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IFU-TLB v01/2019-PT

OBSOLETE

INSTRUCTIONS FOR THE USE OF LAB SCREWS

1. Product description

Lab screws are machined refractory stainless steel (resistant to high temperatures) or medical grade titanium items.
They are longer than the casting cylinder chimney to the wax from wax models entering blocking access.

2. Planned use:

They are used to temporarily anchor casting cylinders to the implant replicate during modelling then to hold the metal coping while building the prosthesis.

3. Directions for use

During prosthesis construction in the prosthetics laboratory they replace the final medical grade titanium clinical screws so that the latter are not damaged in any way.

4. Contraindications

Allergy or hypersensitivity to the chemical components of the materials used.

5. Side effects and interactions

Follow the instructions for use recommended by the manufacturer. Incorrect use of the product may cause injury to the patient.

6. Cleaning, disinfection and sterilisation

In the event of any packaging defect, the product should not be used and the manufacturer notified immediately.

All BIONER lab screws are supplied clean and appropriately packaged. As they are not inserted in the mouth and their use is restricted to the laboratory, they require no sterilisation.

7. Warning

Lab screws should not be reused; they are for single use only as the connections may be damaged, thus altering their technical specifications.

8. Precautions

It is especially important to ensure that each reference is used in combination with dimensionally compatible elements as indicated on the product label.

9. Compatibility


BIONER lab screw family references must be used in conjunction with compatible components. Carefully examine the components available in the company product catalogue for suitability before insertion.

10. Validity

These instructions for use correspond to the current version of references for the lab screw family as indicated at the foot of the document.

11. Storage

BIONER products do not require any special storage conditions even though it is recommended not to expose them to conditions that could affect the product, its packaging or label, such as extreme temperature and humidity, knocks, etc.

ISO 15223 Symbol	Description
 To open	Indicates the direction for opening the blister
	BIONER SA logo
REF	Catalogue reference
LOT AOFAAXXXX	Batch code
 DD/MM/YY	Date of Manufacture
	Refer to the instructions for use
	Do not use if package is damaged
Family	Name of the product family
Mat = Stainless steel/titanium	General indication of the material used
n=1	Indicates the number contained in each blister pack
	Non-sterile
	QR code (reference and batch)

	Do not reuse
	Consult warnings and precautions in instructions for use
	Manufacturer
	EC Marking

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Name BIONER, S.A.